
REF M94 MERRET® TUBES

GB

INTENDED USE

Microgen Bioproducts Merret® Tubes are intended for the staining of reticulocytes in sequestrenated blood samples prior to counting. The Tubes are for professional use by operatives familiar with the preparative techniques for producing blood films and conducting reticulocyte counts.

PRINCIPLE OF THE TEST

Reticulocytes are juvenile red blood cells and contain remnants of ribosomes and ribonucleic acids which were present in larger amounts in the cytoplasm of the nucleated precursors from which they were derived. Ribosomes have the property of reacting with Brilliant Cresyl Blue to form a blue precipitate. The most immature reticulocytes are those with the largest amount of precipitable ribosomal material.

The inner base of each Merret® Tube is coated with a standardised preparation of Brilliant Cresyl Blue which reacts with the reticulocytes in sequestrenated blood after incubation at 37° C for 10 minutes, thus staining the reticulocytes ready for microscopic examination.

The staining procedure is based upon the method of Conn (1940) essentially as originally described by Robertson (1917), who described the staining of reticulocytes using Brilliant Cresyl Blue. The reticulum is stained blue and the rest of the corpuscle is unstained.

CONT

KIT PRESENTATION

96 x Tubes coated with Brilliant Cresyl Blue, protected with an adhesive seal.
Instructions for Use

Materials and Equipment required but not provided

- 37° C Incubator
- Pipette with disposable tips to deliver 250 microlitres or dropper pipette delivering drop size of approximately 60-70 microlitres blood.
- Timer
- Latex gloves, safety glasses and other appropriate protective garments
- Biohazardous waste containers
- Racks for Merret® Tubes
- Vortex Mixer (optional)

WARNINGS AND PRECAUTIONS

Safety

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. All human blood samples should be treated as potentially infectious.
3. Waste reagents and materials should be decontaminated and discarded according to national guidelines. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30minutes. Liquid waste containing acid must be neutralised before treatment.

Procedural

1. The kit should be used in accordance with these Instructions for Use.
2. Do not cross-contaminate samples.

STORAGE AND SHELF LIFE

All unopened components can be used until the date printed on the outer carton label provided they are stored at 15-25°C
Tubes may be removed from the protective seal as required for use. Provided the seal is kept in place remaining tubes may be used up to their expiry date.

Specimens:

Blood should be collected into anti-coagulant and mixed gently. Samples should be stained using Merret® Tubes and counted within 24 hours for best results, although the count will tend to fall as some of the reticulocytes may mature *in vitro*.

TEST PROCEDURE

1. Place 250 microlitres (4 large drops) of sequestrenated blood into the bottom of a Merret® Tube.
2. Mix the blood thoroughly with the stain by holding the top of the tube and vortexing or tapping smartly for at least 10 seconds or until the dried stain has dissolved from the bottom of the tube.
3. Incubate the tube at 37° C for 10minutes
4. Remove a small amount of blood and prepare a blood film
5. Count stained reticulocytes microscopically
6. The reticulum of immature red blood cells (reticulocytes) should stain deep blue.

LIMITATIONS OF USE

1. Results from Reticulocyte counts should always be interpreted in the context of all available clinical and laboratory data.
2. Results from Reticulocyte counts are influenced by the quality of the blood film prepared and visual acuity of the observer. Reticulocyte counts should only be conducted by competent operatives.

PERFORMANCE CHARACTERISTICS

The quality of the stain used in the manufacture of Merret Tubes is monitored to ensure consistency. The standards are set as a UV-VIS peak range for the dye of 624-628nm, minimum dye content 60% and each lot is tested for the staining of reticulocytes.

REPRODUCIBILITY

Intra-batch reproducibility – This has been established by determining the optical densities and staining performance of reconstituted stain over a whole batch of tubes. The co-efficient of variation across the run for stain concentration was <5%

Inter-batch reproducibility – This has been established by determining the optical densities and staining performance of reconstituted stain over three batches of Merret® Tubes. The co-efficient of variation across the three batches for stain concentration was <10%.

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Microgen Bioproducts Ltd