

TRUSTED RESULTS

# Campylobacter Testing

*Offers accurate results, fast*

- Detects more pathogenic species of *Campylobacter*, including *C. jejuni*, *C. coli*, *C. lari*, and *C. upsaliensis*, than any other commercially-available assay
- Provides excellent accuracy and results within minutes without the downsides of culture
- Compatible with multiple specimen types, including fresh, frozen, and samples in transport media



**CAMPYLOBACTER QUIK CHEK™**

**CAMPYLOBACTER CHEK™**

## Campylobacter Test Kits

### Product Order Information

Product Name	Cat. No.	Sample Type	Time to Result	Format
<b>CAMPYLOBACTER QUIK CHEK™</b>	T31025	Fresh, frozen, transport media (Cary Blair or C&S)	30 minutes	Rapid EIA
<b>CAMPYLOBACTER CHEK™</b>	T31096	Fresh, frozen, transport media (Cary Blair or C&S)	60 minutes	Automatable 96-well plate

# CAMPYLOBACTER QUIK CHEK™ Performance

N = 1552	Culture Positive	Culture Negative
<b>CAMPYLOBACTER QUIK CHEK™ Positive</b>	34	13*
<b>CAMPYLOBACTER QUIK CHEK™ Negative</b>	1*	1504

	95% Confidence Limits	
Sensitivity	97.1%	85.5% - 99.9%
Specificity	99.1%	98.5%- 99.5%

\* The 14 discrepant specimens were further characterized by additional testing at TECHLAB. This testing included a US FDA-cleared commercial Microassay well EIA, a US FDA-cleared commercial molecular test, in-house PCR (detecting the 16S rRNA gene of *Campylobacter* spp., and species-specific identification), and bidirectional sequencing. All of the discrepant specimens resolved in favor of the **CAMPYLOBACTER QUIK CHEK™** test.

# CAMPYLOBACTER CHEK™ Performance


N = 1552	Culture Positive	Culture Negative
<b>CAMPYLOBACTER CHEK™ Positive</b>	32	14*
<b>CAMPYLOBACTER CHEK™ Negative</b>	3**	1503

	95% Confidence Limits	
Sensitivity	91.4%	75.8% - 97.8%
Specificity	99.1%	98.5%- 99.5%

The 17 discrepant specimens were further characterized by additional testing at TECHLAB. This testing included a US FDA-cleared commercial Microassay well EIA, a US FDA-cleared commercial molecular test, in-house PCR (detecting the 16S rRNA gene of *Campylobacter* spp., and species-specific identification), and bidirectional sequencing.

\* All of the 14 culture negative and **CAMPYLOBACTER CHEK™** positive resolved in favor of the **CAMPYLOBACTER CHEK™** test.

\*\*One of the three specimens that were culture positive and **CAMPYLOBACTER CHEK™** test negative was confirmed to be negative with all tests.



CAMPYLOBACTER QUIK CHEK™ Package Insert and CAMPYLOBACTER CHEK™ Package Insert.

Distributed by



T: 02 8212 4074 F: 02 9423 6992  
 info@keydiagnostics.com.au  
 www.keydiagnostics.com.au  
 PO Box 1038, Gympie, NSW, 2227

Developed and manufactured by



PROVEN DIAGNOSTIC PERFORMANCE

www.techlab.com

TECHLAB, Inc. is a leading developer and manufacturer of infectious disease diagnostics products worldwide. Since 1989, the company has offered a quality portfolio of diagnostic tests which now includes assays for *C. difficile*, *Campylobacter*, Shiga toxin producing *E. coli*, *H. pylori*, protozoan parasites, and biomarkers of intestinal inflammation.

© 2021 TECHLAB, Inc. All rights reserved.

CAMPYLOBACTER CHEK, CAMPYLOBACTER QUIK CHEK, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc. All other trademarks referenced are trademarks of their respective owners. PN 832021001