


RAPID *C. DIFFICILE* TESTING

C. DIFF QUIK CHEK®

Accurately detect C. difficile antigen (GDH) in fecal specimens

- Guidelines recommend GDH screening in combination with toxin testing
- Algorithm testing provides a more complete diagnostic picture than molecular testing alone
- Molecular assays only indicate if a gene is present, not if disease-causing toxins are being produced
- Sensitivity and Negative Predictive Value (NPV) of GDH are equivalent to toxigenic culture



C. DIFF QUIK CHEK®





Product Order Information

Product Name	Cat. No.	Sample Type	Time to Result	Format
<i>C. DIFF QUIK CHEK</i> ®	30390	Fresh and frozen fecal sample (unpreserved), Samples in Cary Blair or C&S	< 30 minutes	Rapid EIA

C. DIFF QUIK CHEK®


Assay Protocol

- 1 Dilute fecal specimen**
Add 25 µL of specimen to tube containing 500 µL *Diluent* and 1 drop of *Conjugate*. Mix thoroughly.
- 2 Add sample**
Transfer 400 µL of diluted sample-conjugate mixture to *Sample Well*. Incubate at room temperature for 15 minutes.
- 3 Wash Membrane Device**
Add 300 µL *Wash Buffer* to *Reaction Window*. Allow to completely absorb.
- 4 Add Substrate reagent**
Add 2 drops *Substrate* to *Reaction Window*. Incubate at room temperature. Read at 10 minutes after addition of *Substrate*.



Easy-to-read Results

- Results in 30 minutes without instrumentation
- Uses the popular TECHLAB® QUIK CHEK™ format
- Simple and clear interpretation of results
- Refer to the package insert for additional interpretation of results.




Positive



Negative

American Society for Microbiology (ASM). 2010. A Practical Guidance Document for the Laboratory Detection of Toxigenic *Clostridium difficile*.
Crobach MJ, Planche T, Eckert C, Barbut F, Terveer EM, Dekkers OM, et al. European Society of Clinical Microbiology and Infectious Diseases: Update of the diagnostic guidance document for *Clostridium difficile* infection. Clin Microbiol Infect 2016;22:S63e81.
C. DIFF QUIK CHEK Package Insert



Distributed by

keydiagnostics

T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
PO Box 1038, Gymea, NSW, 2227

Developed and manufactured by

TECHLAB®

PROVEN DIAGNOSTIC PERFORMANCE

www.techlab.com

TECHLAB, Inc. is a leading developer and manufacturer of infectious disease diagnostics products worldwide. Since 1989, the company has offered a quality portfolio of diagnostic tests which now includes assays for *C. difficile*, *Campylobacter*, Shiga toxin producing *E. coli*, *H. pylori*, protozoan parasites, and biomarkers of intestinal inflammation.

© 2021 TECHLAB, Inc. All rights reserved.

C. DIFF QUIK CHEK, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc. All other trademarks referenced are trademarks of their respective owners. PN 832021016