

MEANINGFUL RESULTS

***E. HISTOLYTICA II™****Provide straightforward results  
for better patient care*

- Specifically detect *E. histolytica* without cross-reacting with non-pathogenic *E. dispar*
- Avoid inappropriate therapy and unnecessary antibiotic treatment
- Report results in hours, not days


***E. HISTOLYTICA II™*****Product Order Information**

Product Name	Cat. No.	Sample Type	Time to Result	Format
<b><i>E. HISTOLYTICA II™</i></b>	T5017	Fresh and frozen fecal sample (unpreserved)	< 2.5 hours	Microplate ELISA


# E. HISTOLYTICA II™

## Assay Protocol


**1 Dilute fecal specimen**  
Add 400 µL of liquid specimen to 400 µL *Diluent*. Alternatively, for a formed specimen, prepare by using a swab to transfer a sample to 400 µL *Diluent*.




**2 Add Conjugate**  
Add 1 drop (50 µL) of *Conjugate* to each well.




**3 Add samples**  
Add controls or diluted specimen to wells containing *Conjugate*.  
Incubate at room temperature for 2 hours.





**4 Wash Assay Wells**  
Wash wells with diluted *Wash Buffer* 5 times. Slap inverted plate firmly on paper towels between each wash.



**5 Add Substrate reagent**  
Add 2 drops *Substrate* to all test wells. Incubate wells at room temperature for 10 minutes. Tap wells gently to mix during incubation.



**6 Add Stop Solution**  
Add 1 drop *Stop Solution*. Wait 2 minutes before reading. Refer to the package insert for interpretation of results.



E. HISTOLYTICA CHEK Package Insert  
Haque, R. et al. 2003. Current Concepts Amebiasis. N Engl J Med. 348:1565-73.

Distributed by



T: 02 8212 4074 F: 02 9423 6992  
info@keydiagnostics.com.au  
www.keydiagnostics.com.au  
PO Box 1038, Gympie, NSW, 2227

Developed and manufactured by



PROVEN DIAGNOSTIC PERFORMANCE

www.techlab.com

TECHLAB, Inc. is a leading developer and manufacturer of infectious disease diagnostics products worldwide. Since 1989, the company has offered a quality portfolio of diagnostic tests which now includes assays for *C. difficile*, *Campylobacter*, Shiga toxin producing *E. coli*, *H. pylori*, protozoan parasites, and biomarkers of intestinal inflammation.

© 2021 TECHLAB, Inc. All rights reserved.  
E. HISTOLYTICA II, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc. All other trademarks referenced are trademarks of their respective owners. PN 832021008