Petrifilm™

6410/6411/6416/6443

Product Instructions

KEYdiagnostics

T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
PO Box 1038, Gympie, NSW, 2227

-  (EN) Coliform Count Plate
-  (FR) Test pour la numération des Coliformes
-  (DE) Coliforme Zählplatte
-  (IT) Piastra per il conteggio dei coliformi
-  (ES) Placa para recuento de coliformes
-  (NL) Coliform Telplaat
-  (SV) Odlingsplatta för koliforma bakterier
-  (DA) Coliform Tælleplade
-  (NO) For koliforme bakterier
-  (FI) Koliformien kasvatusalusta
-  (PT) Placa para Contagem de Coliformes
-  (EL) Πλακίδιο Καταμέτρησης Κολοβακτηριδίων
-  (PL) Płytko do oznaczania liczby bakterii z grupy coli
-  (RU) Тест-пластина для подсчета колиформных бактерий
-  (TR) Koliform Sayım Plakası
-  (JA) 大腸菌群数測定用プレート
-  (ZH) 大肠菌群测试片
-  (TH) แผ่นอาหารเลี้ยงเชื้อสำหรับนับจำนวนโคลิฟอร์ม
-  (KO) 대장균군 측정용 플레이트

Coliform Count

Product Instructions

Coliform Count Plate

Product Description and Intended Use

The 3M™ Petrifilm™ Coliform Count (CC) Plate is a sample-ready-culture medium system which contains modified Violet Red Bile (VRB) nutrients, a cold-water-soluble gelling agent, and a tetrazolium indicator that facilitates colony enumeration. 3M Petrifilm CC Plates are used for the enumeration of coliforms in the food and beverage industries. 3M Petrifilm CC Plate components are decontaminated though not sterilized. 3M Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. 3M Petrifilm CC Plate has not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

The United States Food and Drug Administration (FDA) Bacteriological Analytical Manual (BAM)⁷ define coliforms as Gram-negative rods, which produce acid and gas from lactose fermentation. Coliform colonies growing on the 3M Petrifilm CC Plate produce acid, which causes the pH indicator to deepen the gel color; gas trapped around red colonies indicates coliforms.

ISO defines coliforms by their ability to grow in method-specific, selective media. ISO method 4832¹, enumerating coliforms by colony-count technique, defines coliforms as acid producers on VRB with lactose (VRBL) agar. On 3M Petrifilm CC Plates these acid-producing coliforms are indicated by red colonies with or without gas production. ISO method 4831², enumerating coliforms by the most probable number (MPN) method, defines coliforms by their ability to grow and produce gas from lactose in a selective broth. On 3M Petrifilm CC Plates these coliforms are indicated by red colonies associated with gas.

AFNOR Certification has certified 3M Petrifilm CC Plate in comparison to ISO method 4831² and ISO method 4832¹ for enumeration of total coliforms. AFNOR Certification has also certified 3M Petrifilm CC Plate in comparison to NF V08-060³, for enumeration of thermotolerant coliforms.

Safety

The user should read, understand, and follow all safety information in the instructions for the 3M Petrifilm CC Plate. Retain the safety instructions for future reference.

⚠ WARNING: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instruction contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform 3M Petrifilm CC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices⁴, ISO 7218⁵, or ISO 17025⁶.

To reduce the risks associated with misinterpretation of results:

- 3M has not documented 3M Petrifilm CC Plates for use in industries other than food and beverage. For example, 3M has not documented 3M Petrifilm CC Plates for testing water, pharmaceuticals, or cosmetics.
- Do not use the 3M Petrifilm CC Plates in the diagnosis of conditions in humans or animals.
- 3M Petrifilm CC Plates do not differentiate any one coliform strain from another.
- Foods with high sugar content may increase the potential for gas production from non-coliform *Enterobacteriaceae*.

Consult the Safety Data Sheet for additional information.

For information on documentation of product performance, visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its

authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

Limitation of 3M Liability

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Storage

Store unopened 3M Petrifilm CC Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened 3M Petrifilm CC Plate pouches to come to room temperature before opening. Return unused 3M Petrifilm CC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of 3M Petrifilm CC Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place 3M Petrifilm CC Plates in a sealable container. To remove frozen 3M Petrifilm CC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container. 3M Petrifilm CC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the 3M Petrifilm CC Plates to moisture which can damage the plates.

Do not use 3M Petrifilm CC Plates that show discoloration. Expiration date and lot number are noted on each package of 3M Petrifilm CC Plates. The lot number is also noted on individual 3M Petrifilm CC Plates.

⚠ Disposal

After use, 3M Petrifilm CC Plates may contain microorganisms that may be a potential biohazard.

Follow current industry standards for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffered dilution water⁷, 0.1% peptone water⁸, peptone salt diluents^{8,9}, Buffered Peptone Water^{8,9}, saline solution (0.85-0.90%), bisulfite-free letheen broth or distilled water. See section "Specific Instructions for Validated Methods" for specific requirements.

Do not use diluents containing citrate, bisulfite or thiosulfate with 3M Petrifilm CC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C (104-113°F).

2. Blend or homogenize sample.
3. For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6.6 - 7.2. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

Plating

1. Place the 3M Petrifilm CC Plate on a flat, level surface.
2. Lift the top film and with the pipette perpendicular to the inoculation area dispense 1 mL of sample suspension onto the center of bottom film.

3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the 3M™ Petrifilm™ Spreader with the flat side down on the center of the 3M Petrifilm CC Plate. Press gently on the center of the 3M Petrifilm Spreader to distribute the sample evenly. Spread the inoculum over the entire 3M Petrifilm CC Plate growth area before the gel is formed. Do not slide the 3M Petrifilm Spreader across the film.
5. Remove the 3M Petrifilm Spreader and leave the 3M Petrifilm CC Plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate 3M Petrifilm CC Plates in a horizontal position with the clear side up in stacks of no more than 20 plates. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the “**Specific Instructions for Validated Methods**” section.

Interpretation

1. 3M Petrifilm CC Plates can be counted using a standard colony counter or other illuminated magnifier. Gas produced by coliform may disrupt the colony so that the colony “outlines” the bubble. This should be counted as a single coliform. Do not count colonies on the foam dam since they are removed from the selective influence of the medium. Do not count artifact bubbles that may be present.

The interpretation of coliform colonies on the 3M Petrifilm CC Plate varies by reference method. For example:

AOAC® Official MethodsSM:

Coliform colonies are red and closely associated (within one colony diameter) with entrapped gas. Colonies not associated with gas (a distance greater than one colony diameter between colony and gas bubble) are not counted as coliforms.

OR NF Validation certified methods:

- As compared to ISO method 4831² (MPN method), coliform colonies are red and closely associated (within one colony diameter) with entrapped gas.
- As compared to ISO method 4832¹ (total coliforms) and NF V08-060³ (fecal coliforms) both based on VRBL methods, count all red colonies with or without gas.

2. The circular growth area is approximately 20 cm². Estimates can be made on 3M Petrifilm CC Plates containing greater than 150 colonies by counting the number of colonies in one or more representative squares and determining the average number per square. Multiply the average number by 20 to determine the estimated count per 3M Petrifilm CC Plate.
3. When colonies are present in large numbers, 3M Petrifilm CC Plates may have a deepening of the gel color and either or both of the following characteristics: many small, indistinct colonies and/or many gas bubbles. High concentrations of coliforms will cause the growth area to turn dark red. When this occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
4. Where necessary, colonies may be isolated for further identification. Lift the top film using proper testing technique and pick the colony from the gel. Test using standard procedures.
5. If the 3M Petrifilm CC Plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C (5°F) for no longer than one week.

For further information refer to the appropriate “3M™ Petrifilm™ Coliform Count Plate Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

Specific Instructions for Validated Methods

AOAC® Official MethodsSM (986.33 Bacteria and Coliform Counts in Milk, Dry Rehydratable Film Methods and 989.10 Bacterial and Coliforms Counts in Dairy Products, Dry Rehydratable Film Methods)

Incubate 3M Petrifilm CC Plates 24 hours ± 2 hours at 32°C ± 1°C.

AOAC® Official MethodsSM (991.14 Coliform and *Escherichia coli* Counts in Foods, Dry Rehydratable Film Methods)

Incubate 3M Petrifilm CC Plates 24 hours ± 2 hours at 35°C ± 1°C.

NF Validation by AFNOR Certification:

NF Validation certified method in compliance with ISO 16140¹⁰ in comparison to ISO 4832¹ (3M-01/2-09/89 A)

Scope of the validation: For testing all human food products (except raw shellfish), pet food and environmental samples.

Use the following details when implementing the above Instructions for Use:

Sample preparation:

Use only ISO listed diluents.

Incubation:

Incubate 3M Petrifilm CC Plates 24 hours ± 2 hours at 30°C ± 1°C or 37°C ± 1°C.

Interpretation:

Calculate the number of microorganisms present in the test sample according to ISO 7218⁵ for one plate per dilution. Estimates are outside of the scope of the NF Validation certification (cf interpretation part paragraph 2).

**NF Validation certified method in compliance with ISO 16140¹⁰ in comparison to ISO 4831² (3M-01/2-09/89 B)**

Scope of the validation: For testing all human food products (except raw shellfish).

Use the following details when implementing the above Instructions for Use:

Sample preparation:

Use only ISO listed diluents.

Incubation:

Incubate 3M Petrifilm CC Plates 24 hours ± 2 hours at 30°C ± 1°C or 37°C ± 1°C.

Interpretation:

Calculate the number of microorganisms present in the test sample according to ISO 7218⁵ for one plate per dilution. Estimates are outside of the scope of the NF Validation certification (cf interpretation part paragraph 2).

NF Validation certified method in compliance with ISO 16140¹⁰ in comparison to NF V08-060³ (3M-01/2-09/89 C)

Scope of the validation: For testing all human food products

Use the following details when implementing the above Instructions for Use:

Sample preparation:

Use only ISO listed diluents

Incubation:

Incubate 3M Petrifilm CC Plates 24 hours ± 2 hours at 44°C ± 1°C.

Interpretation:

Calculate the number of microorganisms present in the test sample according to ISO 7218⁵ for one plate per dilution. Estimates are outside of the scope of the NF Validation certification (cf interpretation part paragraph 2).

**3M 01/2 – 09/89A, 01/2 – 09/89B, 01/2 – 09/89C
ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS**

<http://nf-validation.afnor.org/en>

For more information about end of validity, please refer to NF VALIDATION certificate available on the website mentioned above.

References

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms –Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

Refer to the current versions of the standard methods listed above.

EXPLANATION OF SYMBOLS

www.3M.com/foodsafety/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL
Official Methods is a service mark of AOAC INTERNATIONAL

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Europe and MEA

3M Deutschland GmbH
Carl-Shurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M United Kingdom PLC

Morley Street, Loughborough
Leicestershire
LE11 1EP
United Kingdom
+(44) 1509 611 611

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3M Health Care

2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

© 2015, 3M. All rights reserved.
3M and Petrifilm are trademarks of 3M. Used under license in Canada.
34-8717-3405-8

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Références

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
 2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
 3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
 4. U.S. Food and Drug Administration Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
 5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
 6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
 7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
 8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
 9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
 10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.
- Se reporter aux versions en cours de validité des méthodes normalisées citées plus haut.

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3M 01/2 – 09/89A, 01/2 – 09/89B, 01/2 – 09/89C
ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

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Literaturangaben

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM gefunden unter: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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Riferimenti

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
 2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
 3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
 4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
 5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
 6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
 7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
 8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
 9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
 10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.
- Fare riferimento alle versioni attuali dei metodi standard elencati in precedenza.

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Referencias

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Manual de análisis bacteriológico (BAM), 8va. Edición, Revisión A, 1998. Encuentre el índice de reactivos para (BAM) en: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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3M 01/2 – 09/89A, 01/2 – 09/89B, 01/2 – 09/89C
ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS
<http://nf-validation.afnor.org/en>

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Referenties

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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Referenser

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagensindex för BAM återfinns på: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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Referanser

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagentindeks for BAM ble funnet på: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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Referências

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44 °C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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3M 01/2 – 09/89A, 01/2 – 09/89B, 01/2 – 09/89C ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

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1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM που βρίσκεται στο: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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Źródła

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
 2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
 3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
 4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
 5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
 6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
 7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Lista odczynników dla BAM: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
 8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
 9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
 10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.
- Należy odnieść się do aktualnych wersji wymienionych powyżej metod standardowych.

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Referanslar

1. ISO 4832. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831. Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms –Most probable number technique.
3. AFNOR NF V08-060 General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM şu adreste mevcuttur: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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1. ISO 4832 Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
2. ISO 4831 Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms –Most probable number technique.
3. AFNOR NF V08-060. General guidance for the enumeration of fecal coliforms – Colony count technique (VRBL) at 44°C – Routine method.
4. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
5. ISO 7218 Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. ISO/IEC 17025 General requirements for the competence of testing and calibration laboratories.
7. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>
8. ISO 6887-5. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 5: Specific rules for the preparation of milk and milk products.
9. ISO 6887-1 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination, Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
10. ISO 16140 Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods.

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3M Österreich GmbH
3M Health Care Limited