**Compact Dry LS** medium for *Listeria spec* / Medium zum Nachweis von *Listerien* / milieu pour *listeria* / medio para *listeria* / medio per *listeria* / meio para *listeria*

40 plates/Platten/plaques/placas/piastre/placas

ID-No. 1 000 899

240 plates /Platten/plaques/placas/piastre/placas

ID-No. 1 010 013

920 plates/Platten/plaques/placas/piastre/placas

ID-No. 1 002 881

### English

### Deutsch

### Français

*Compact Dry LS is a ready to use, chromogenic plate for detection of Listeria spec.*

#### Sample pretreatment

**Viable count in water or liquid foodstuff**

Drop 1 ml of specimen (dilute if necessary) on the middle of the Compact Dry plate.

#### Viable count in solid foodstuff

Add Buffered Peptone Water (BPW) buffer to the sample and homogenize by stomacher®. Drop 1 ml of specimen (dilute if necessary) on the middle of the dry sheet of the Compact Dry plate.

#### Viable count in swab test specimen

Use the swab to wipe the surface, put into the device with wiping solution. Drop 1 ml of wiping solution (dilute if necessary) on the middle of the Compact Dry plate. It is recommended to use "Swab for Compact Dry" offered by HyServe Id-No. 1 002 952/3 (40/240 pieces).

#### Test instructions

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator.
5. After incubation count the number of colored colonies underneath the plate. White paper placed under the plate helps to count the colonies.

**Incubation time** 24 ± 2 hours  
plus a further 24 ± 2 hours

**Incubation temperature** 35 or 37 ± 1 °C

#### Interpretation of the results

*Listeria* species form light blue/blue colonies of 1-2 mm in diameter.

#### Storage and shelf life

Keep at room temperature (+1 to +30 °C). Total shelf life 12 months after manufacturing.

#### Notes

- It is known that *Listeria ivanovii* tends to grow slowly and *Listeria seeligeri* tends to be inhibited in this plate.
- A combination of selective agents in the medium inhibits growth of other bacteria.
- High concentrations on plates (> 300 CFU) will cause the entire growth area to become blue/green. In this case dilute the sample.
- After use please follow the current disposal regulations.
- For effective recovery of *Listeria*, it is recommended that samples should be inoculated after an 1 hour resuscitation step at 20 °C by peptone buffered water (BPW) (Hyserve Code 1000889).
- This resuscitation step does not affect the bacterial number and is according to ISO 11290-2, a horizontal method for detection and enumeration of *Listeria monocytogenes*.
- The growth area is 20 cm<sup>2</sup>. The back of the plate has a grid carved of 1 cm x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates are produced at an ISO 9001 certified site.

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