**Compact Dry VP medium for *Vibrio parahaemolyticus***

40 plates/Platten/plaques/placas/lastre/placas	ID-No. 1 000 900
240 plates/Platten/plaques/placas/lastre/placas	ID-No. 1 001 014
500 plates/Platten/plaques/placas/lastre/placas	ID-No. 1 401 014
1400 plates/Platten/plaques/placas/lastre/placas	ID-No. 1 502 886

**English****Deutsch****Français**

*Compact Dry VP is a ready to use, chromogenic plate for detection of Vibrio parahaemolyticus*

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**Sample pretreatment****Viable count in water or liquid foodstuff**

Drop 1 ml of specimen (dilute if necessary) on the middle of the Compact Dry plate.

**Viable count in solid foodstuff**

Add buffer solution to the sample and homogenize by stomacher®. Drop 1 ml of specimen (dilute if necessary) on the middle of the dry sheet of the Compact Dry plate.

**Viable count in swab test specimen**

Use the swab to wipe the surface, put into the device with wiping solution. Drop 1 ml of wiping solution (dilute if necessary) on the middle of the Compact Dry plate. It is recommended to use "Swab for Compact Dry" offered by HyServe Id-No. 1 002 952/3 (40/240 pieces)

**Test instructions**

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator.
5. After incubation count the number of colored colonies underneath the plate. White paper placed under the plate helps to count the colonies.

**Incubation time** 20 - 24 hours  
**Incubation temperature** 35 - 37 °C

Please use the incubation time/temperature according to the national food analysis recommended for vibrios viable count.

**Interpretation of the results**

*Vibrio parahaemolyticus* develops blue/green or blue colonies. Whereas other vibrios develop white (colorless) colonies.

**Storage and shelf life**

Keep at room temperature (+ 5 - 30 °C). Total shelf life 16 months after manufacturing.

**Notes**

- After use please follow the current disposal regulations.
- High concentrations on plates (> 300 CFU) will cause the entire growth area to become blue/green or blue. In this case dilute the sample.
- The growth area is 20 cm<sup>2</sup>. The back of the plate has a grid carved of 1 cm x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates are produced at an ISO 9001 certified site.

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