

## Performance Evaluation of InSite *Listeria* (IL050/IL100) and New InSite *Listeria* (ILC050/ILC100)

### Introduction

InSite *Listeria* is a screening test for *Listeria* spp., intended for the use on food contact surfaces and food processing equipment after cleaning to detect the presence of *Listeria* species. Insite *Listeria* is a self-contained swab device containing a selective growth medium to inhibit most non-*Listeria* microorganisms and provide nutrients to support growth of *Listeria*. Indicator compounds turn broth from yellow to black utilizing B-glucosidase enzyme production produced by *Listeria* species. A darkening of the media to grey/black after 24-48 hours at 37°C indicates a presumptive positive result for *Listeria* spp.

A new InSite *Listeria* Species device has been developed to improve selectivity reducing the proportion of presumptive positives from non-*Listeria* that can survive and grow in the current InSite *Listeria*. This study was designed to evaluate the performance of detecting *Listeria* species as well as selectivity against non-*Listeria* species.

### Equipment, Supplies and Reagents

- InSite *Listeria*- IL050
- New InSite *Listeria* – ILC050
- Hygiena™ Digital Dry Block Incubator for maintaining temperatures at 37 ± 2°C – INCUBATOR2
- Hygiena Swab Tube Block for Model No: INCUBATOR2 – IB001
- Brain heart infusion broth (BHI)
- Buffered Peptone Water (BPW)

Table 1. Inclusivity Panel for <i>Listeria</i> species	
Sample No.	Genus species
1	<i>L. monocytogenes</i> ATCC 19115
2	<i>L. monocytogenes</i> ATCC 19118
3	<i>L. monocytogenes</i> ATCC 7644
4	<i>L. monocytogenes</i> ATCC BAA-751
5	<i>L. monocytogenes</i> ATCC 13932
6	<i>L. innocua</i> NCTC 11288
7	<i>L. innocua</i> ATCC 33090
8	<i>L. welshimeri</i> ATCC 35897
9	<i>L. seeligeri</i> ATCC 35967
10	<i>L. ivanovii</i> ATCC 19119

**Table 2. Exclusivity Panel for Non-*Listeria* species <sup>a</sup>**

Sample No.	Genus species	Sample No.	Genus species
1	<i>Klebsiella pneumoniae</i>	13	Unknown ID #1
2	<i>Bacillus licheniformis</i>	14	<i>Bacillus licheniformis</i>
3	<i>Pseudomonas fragi</i>	15	<i>Proteus vulgaris</i>
4	<i>Bacillus licheniformis</i>	16	Unknown ID #2
5	<i>Enterococcus faecalis</i>	17	<i>Bacillus pumilus</i>
6	<i>Carnobacterium maltaromaticum</i>	18	<i>Proteus vulgaris</i>
7	<i>Carnobacterium maltaromaticum</i>	19	<i>Bacillus pumilus</i>
8	<i>Bacillus pumilus</i>	20	<i>Enterococcus faecalis</i>
9	<i>Micrococcus luteus</i>	21	<i>Streptomyces cinnamoneus</i>
10	<i>Streptomyces cinnamoneus</i>	22	<i>Carnobacterium maltaromaticum</i>
11	<i>Klebsiella pneumoniae</i>	23	<i>Bacillus pumilus</i>
12	<i>Bacillus licheniformis</i>	24	<i>Bacillus pumilus</i>

<sup>a</sup> Confirmed Non-*Listeria* Organisms Isolated from Industrial Environmental samples. ID using Hygiena RiboPrinter® System

## Sample Preparation and Enrichment

**Inclusivity Panel** – Inclusivity cultures were enriched by adding a colony from a TSA plate to BHI broth for 24 hours producing an overnight culture. Each culture was then serially diluted 1:10 to produce a concentration between 10<sup>3</sup>-10<sup>4</sup>CFU/mL. 100 µL were spread onto TSA plates to perform a purity check and confirm final CFU levels. For each suspension 10µL was pipetted onto the top of duplicate swabs producing a concentration of 10<sup>1</sup>-10<sup>2</sup> CFU per swab. Devices were activated and incubated at 37°C for 48 hours.

**Exclusivity Panel** – Exclusivity cultures were enriched by adding a colony from a TSA plate to BHI broth for 24 hours producing an overnight culture. Each Culture was then serially diluted 1:10 to -1, -2, and -3. 100 µL were spread onto TSA plates to perform a purity check and confirm final CFU levels. For each suspension neat to -3, 10µL was pipetted onto the top of swab buds. Devices were activated and incubated at 37°C for 48 hours.

## Method

For each device, the media was visually inspected for color change from yellow/amber to black/grey after 48hrs. Results were recorded as presumptive positive or negative for all samples. All samples were confirmed according to BAM chapter 10, Detection of *Listeria* monocytogenes in Foods and Environmental Samples (2).

## Results and Discussion

InSite *Listeria* and New InSite *Listeria* devices demonstrated a presumptive positive color change for all *Listeria* species tested (Table 3). This is in line with the historical AOAC RI certification for InSite media and the AOAC RI certification of New InSite *Listeria* (1,3).

The exclusivity results of InSite *Listeria* and New InSite *Listeria* were compared using probability of detection (POD) and difference in POD (dPOD) (Table 4). Of the non-*Listeria* samples tested, both InSite *Listeria* and New InSite *Listeria* detected non-*Listeria* species at a CFU per swab greater than 1000 CFU per swab. As the concentration of non-*Listeria* species increased the presumptive positive POD increases for both devices. The statistical analysis demonstrated a significant difference between the two tests when each test contained greater than  $10^3$  CFU per swab.

Table 3. Inclusivity Results IL50 vs ILC050 <sup>bc</sup>			
Sample No.	Genus species	InSite <i>Listeria</i> (IL50) Result <sup>a</sup>	New InSite <i>Listeria</i> (ILC050) Result <sup>a</sup>
1	<i>L. monocytogenes</i> ATCC 19115	+	+
2	<i>L. monocytogenes</i> ATCC 19118	+	+
3	<i>L. monocytogenes</i> ATCC 7644	+	+
4	<i>L. monocytogenes</i> ATCC BAA-751	+	+
5	<i>L. monocytogenes</i> ATCC 13932	+	+
6	<i>L. innocua</i> NCTC 11288	+	+
7	<i>L. innocua</i> ATCC 33090	+	+
8	<i>L. welshimeri</i> ATCC 35897	+	+
9	<i>L. seeligeri</i> ATCC 35967	+	+
10	<i>L. ivanovii</i> ATCC 19119	+	+

<sup>a</sup> Chromogenic color change result

<sup>b</sup>  $10^1$ - $10^2$  CFU per swab

<sup>c</sup> N=2 for each sample

**Table 4. Exclusivity Results IL050 vs ILC050**

Sample Type	CFU/Test Portion	N	InSite <i>Listeria</i> (IL050)			New InSite <i>Listeria</i> (ILC050)			dPOD	95% CI
			X	POD <sub>IL</sub>	95% CI	X	POD <sub>ILc</sub>	95% CI		
Non- <i>Listeria</i> Species	3.8x10 <sup>7</sup>	4	4	1.00	0.51, 1.00	1	0.25	0.00, 0.70	0.75	0.51, 0.30
	4.5x10 <sup>6</sup>	11	9	0.82	0.52, 0.95	4	0.36	0.15, 0.65	0.45	0.37, 0.30
	4.0x10 <sup>5</sup>	19	13	0.68	0.46, 0.85	4	0.21	0.09, 0.43	0.47	0.38, 0.41
	4.0x10 <sup>4</sup>	24	12	0.50	0.31, 0.69	2	0.08	0.02, 0.26	0.42	0.29, 0.43
	4.1x10 <sup>3</sup>	20	6	0.30	0.15, 0.52	2	0.10	0.03, 0.30	0.20	0.12, 0.22
	3.6x10 <sup>2</sup>	13	0	0.00	0.00, 0.23	0	0.00	0.00, 0.23	0.00	0.00, 0.00
	2.7x10 <sup>1</sup>	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	0.00, 0.00

CFU/Test portion = Average CFU per swab

N = Number of test portions

X = Number of positive test portions

POD<sub>IL</sub> = IL50 method positive results divided by the total number of test portions

POD<sub>ILc</sub> = ILC050 method positive results divided by the total number of test portions

dPOD = Difference between IL50 and ILC050 POD values

95% CI = If the confidence interval of dPOD does not contain zero, then the difference is statistically significant at the 5% level

## Conclusions

The results of this study demonstrate that the overall performance of the new InSite *Listeria* (ILC050) outweighs the performance of InSite *Listeria* (IL050). Both devices were able to detect a variety of different *Listeria* species within 48 hours of enrichment. New InSite *Listeria* was able to better suppress the growth of a wide concentration of Non-*Listeria* per swab resulting in fewer presumptive positives. This increased selectivity reduced the number of presumptive positives from non-*Listeria* that would ultimately confirm negative when further confirmation procedures are followed.

## References

1. Calderon, D., Familiari, N., and Meighan, P., InSite *L. mono* Glo for Detection of *Listeria* species and *Listeria monocytogenes* from Environmental Surfaces, AOAC® *Performance Tested*<sup>SM</sup> certification number 121902
2. U.S. Food and Drug Administration (2017) *FDA Bacteriological Analytical Manual*, Chapter 10, Detection of *Listeria monocytogenes* in Foods and Environmental Samples, and Enumeration of *Listeria monocytogenes* in Foods
3. Yurttas, H.C., Maher, J., Danter, W., Bargoo, L., Brown, M., Olstein, A., and Feirtag, J., Evaluation of the PDX-LIB, AOAC® *Performance Tested*<sup>SM</sup> certification number 040501

**key**diagnostics

T: 02 8212 4074 F: 02 9423 6992  
info@keydiagnostics.com.au  
www.keydiagnostics.com.au  
PO Box 1038, Gymea, NSW, 2227