


WaterShot™ | Universal Water ATP Test




Free

Total

WaterShot is a universal ATP water test designed to work with 3M™ or BioTrace™ luminometers. Substituting WaterShot helps to reduce costs and achieve more accurate results. The device is available in two formats: Free and Total. WaterShot Free measures dissolved ATP that is free in solution (non-microbial ATP). WaterShot Total measures both free ATP in solution and microbial ATP (non-microbial and microbial ATP). The difference between Total and Free provides a measure of microbial contamination in the sample. Both tests feature a specifically designed dipper tip that collects 100 µl of water, ensuring consistent sample collection. The test device's unique liquid-stable chemistry, combined with simple design, offers unrivaled accuracy and reproducibility. The devices are often used for the management of CIP systems, water treatment systems such as cooling towers and closed water systems.



Find support documents, instructional videos, and more at www.hygiena.com



Benefits:

- Costs up to 50% less than CleanTrace™ water ATP devices
- Patented Snap-Valve™ design enables accurate results with less variation
- Unique liquid-stable reagent provides accuracy and reproducibility
- Reagent is tolerant to temperature abuse and sanitizer effects
- Honeycomb dipper collects consistent 100 µl sample
- Compatible with 3M™ and BioTrace™ instruments
- 100% recyclable plastic
- 12 month shelf life at refrigerated temperatures (2-8°C)
- 4-week shelf life at room temperatures (21-25°C)

Sensitivity:

- Extremely sensitive (2 fmol ATP)
- Detects down to 10⁴ CFU

| Catalog No. | Description | Quantity |
|----------------|-----------------|----------|
| HG-W-SPXL1333 | WaterShot Total | 100 |
| HG-WF-SPXL1333 | WaterShot Free | 100 |

Trade Up Program

Hygiena's trade-up program allows you to trade in a 3M™ luminometer for a brand new Hygiena luminometer. Our instruments are multi-platform systems, allowing you to collect, analyze, and report data from multiple quality indicators. Unlike other monitoring system manufacturers, there is no scheduled yearly maintenance and calibration necessary for the system. Hygiena's trade-up program can help you save up to 50% on testing costs and thousands of dollars per year.




Visit www.hygiena.com to learn more about Trade Up and see how much you could save!



Related Products:


SnapShot™

SnapShot universal ATP surface tests are designed to be compatible with 3M™ or BioTrace™ luminometers. SnapShot utilizes Hygiena's unique liquid-stable reagent and patented Snap-Valve™ design to deliver more reliable, consistent, and sensitive results than CleanTrace™ test devices. Performance data is extracted from a third part comparative study. Contact us for a full copy.



*All performance data is extracted from a third party comparative study. Contact us for a complete copy.

Repeatability Comparison: Variability of manufacturer devices



| Catalog No. | Description | Quantity |
|-------------|--|----------|
| HG-SPXL1333 | SnapShot Universal ATP for Surfaces-3M | 100 |