

GlutenTox Sticks

QUICK AND RELIABLE

Gluten detection kit for food manufacturers Biomedal GlutenTox Sticks is a rapid gluten detection kit for food and drinks but also for working surfaces. GlutenTox Sticks contains the G12 anti-gliadin antibody that specifically recognizes the most immunogenic fraction of gluten that triggers Celiac Disease.

GlutenTox Sticks advantages

- » Detection limit 3ppm (< 20ppm is considered "gluten free")
- » Results in 1 hour for food samples and 15 minutes for working surfaces
- » Suitable for most hydrolysed food
- » Specific to the toxic fraction of gluten from wheat, barley, rye and oat
- » No cross-reactivity with soy, rice or corn
- » Unique method to detect oat toxicity



GlutenTox Sticks is a practical tool for the control of gluten content in your products.

Each kit contains 25 tests to perform gluten analysis.

Industrial application

- » Semi-quantification of gluten in food
- » Gluten-free food quality control
- » Trace gluten contamination in food and working surfaces
- » Monitoring of gluten detoxification
- » Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



WORKING ON FOOD SAFETY

keydiagnostics

T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
PO Box 1038, Gymea, NSW, 2227



info@biomedal.com
www.biomedal.com

