

# Actero™ Elite **Listeria** Enrichment Media advantages.

FoodChek™'s patented Actero™ Elite Listeria Enrichment Media is a selective bacterial culture medium specifically optimized for the recovery of *Listeria monocytogenes* and *Listeria* spp. and from environmental and food samples in a single-step enrichment. We continue to strive to be the food industry's single, trusted provider for pathogen and other microorganism testing supplies.

Actero™ Elite Listeria Enrichment Media is compatible with the BAX® System Real-Time PCR Assays for *Listeria* and represents the result of true food safety innovation.



**Single-step enrichment**

As short as  
**20 HOURS**

**"time-to-results"**



**Accurate**

-  Reduce sample handling
-  Use less enrichment media\*
-  Improve workflow
-  Higher throughput
-  Release food products sooner

\* Depending on matrix

## Single-Step Enrichment of Environmental and Food Samples with Actero™ Elite Listeria Enrichment Media



# Actero® Elite **Listeria** Enrichment Media product list.

Micro-organism	FoodChek Cat. #	HYGIENA Cat. #	Product Name	Format	Description	Availability		
						United States	Canada	Europe
<b>Listeria</b>	FCM-011	<b>MED2017</b>	ACTERO™ ELITE LISTERIA ENRICHMENT MEDIA	500G	Specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples. 500G <b>dehydrated</b> media (yields 9.29L rehydrated media).			
	FCM-022	<b>MED2018</b>	ACTERO™ ELITE LISTERIA ENRICHMENT MEDIA	2KG	A selective medium specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples. Four (4) 2KG of <b>dehydrated</b> Actero™ Elite Listeria Enrichment Media.			
	FCM-023	<b>MED2020</b>	ACTERO™ ELITE LISTERIA ENRICHMENT MEDIA	10KG	A selective medium specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples. 10KG of <b>dehydrated</b> Actero™ Elite Listeria Enrichment Media.			
	FCM-045	<b>MED2021</b>	ACTERO™ ELITE LISTERIA ENRICHMENT – MEDIABOX™	5L	A selective medium specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples – in <b>liquid</b> format.			
	FCM-046	<b>MED2022</b>		10L	A selective medium specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples – in <b>liquid</b> format.			

2018-08-08F

🕒 Single-step enrichment

🕒 Fastest “time-to-results”

🕒 Targeted accuracy



**key**diagnostics

T: 02 8212 4074 F: 02 9423 6992  
info@keydiagnostics.com.au  
www.keydiagnostics.com.au  
PO Box 1038, Gymea, NSW, 2227

