

ATP Positive Control Kit



Find support documents, instructional videos, and more at www.hygiena.com

The ATP Positive Control Kit validates the efficacy and stability of Hygiena ATP test devices. The kit comes with 25 sealed glass vials containing lyophilized adenosine triphosphate (ATP) and sugars to provide predictable results when sampled by test devices.

Benefits:

- Validates performance of Hygiena test devices
- Ready to use out of the box
- No lab equipment needed
- No reconstitutions or complicated sample preparation
- Meets GMP and other Quality program requirements
- 24 month shelf life at refrigerated temperature (2-8°C)

Compatible products:

- UltraSnap ATP Surface Test
- SuperSnap High Sensitivity ATP Surface Test
- SnapShot Universal ATP Surface Test
- SpotCheck Plus Food Safety Residue Test

| Catalog No. | Description | Quantity |
|-------------|---|----------|
| CK25 | ATP Positive Control Kit | 1 kit |
| US2020 | UltraSnap Surface ATP Test | 100 |
| SUS3000 | SuperSnap High Sensitivity ATP Surface Test | 100 |
| See catalog | SnapShot Universal ATP Test | 100 |
| GL100 | SpotCheck Plus | 100 |



Calibration Control Kit

SystemSURE Plus and EnSURE luminometers run a self-calibration check at start up. However, when additional verification is required, the Calibration Control Kit confirms the luminometer is within specifications and operating correctly. The Calibration Control Kit consists of reusable positive and negative devices. In-house calibration checks confirm systems are working properly and demonstrate due diligence that the monitoring program is in control.



Find support documents, instructional videos, and more at www.hygiena.com

Benefits:

- Compatible with EnSURE or SystemSURE Plus luminometers
- Verifies instrument calibration in-house
- Demonstrates due diligence
- Reusable for life of product
- 5 year shelf life



T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
 PO Box 1038, Gymea, NSW, 2227

| Catalog No. | Description | Quantity |
|-------------|-------------------------|----------|
| PCD4000 | Calibration Control Kit | 1 |
| CK25 | Positive Control Kit | 25 vials |